



Determination of Total Carbohydrates in Algal Biomass

Laboratory Analytical Procedure (LAP)

Issue Date: December 2, 2013

Revision Date: December 29, 2015

Stefanie Van Wychen and Lieve M.L. Laurens

**NREL is a national laboratory of the U.S. Department of Energy
Office of Energy Efficiency & Renewable Energy
Operated by the Alliance for Sustainable Energy, LLC**

This report is available at no cost from the National Renewable Energy Laboratory (NREL) at www.nrel.gov/publications.

Technical Report
NREL/TP-5100-60957
Revised December 2015

Contract No. DE-AC36-08GO28308



Determination of Total Carbohydrates in Algal Biomass

Laboratory Analytical Procedure (LAP)

Issue Date: December 2, 2013

Revision Date: December 29, 2015

Stefanie Van Wychen and Lieve M.L. Laurens

Prepared under Task Nos. BB07.9H10 and BB07.9G10

**NREL is a national laboratory of the U.S. Department of Energy
Office of Energy Efficiency & Renewable Energy
Operated by the Alliance for Sustainable Energy, LLC**

This report is available at no cost from the National Renewable Energy
Laboratory (NREL) at www.nrel.gov/publications.

National Renewable Energy Laboratory
15013 Denver West Parkway
Golden, CO 80401
303-275-3000 • www.nrel.gov

Technical Report
NREL/TP-5100-60957
Revised December 2015

Contract No. DE-AC36-08GO28308

NOTICE

This report was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or any agency thereof.

This report is available at no cost from the National Renewable Energy Laboratory (NREL) at www.nrel.gov/publications.

Available electronically at SciTech Connect <http://www.osti.gov/scitech>

Available for a processing fee to U.S. Department of Energy and its contractors, in paper, from:

U.S. Department of Energy
Office of Scientific and Technical Information
P.O. Box 62
Oak Ridge, TN 37831-0062
OSTI <http://www.osti.gov>
Phone: 865.576.8401
Fax: 865.576.5728
Email: reports@osti.gov

Available for sale to the public, in paper, from:

U.S. Department of Commerce
National Technical Information Service
5301 Shawnee Road
Alexandria, VA 22312
NTIS <http://www.ntis.gov>
Phone: 800.553.6847 or 703.605.6000
Fax: 703.605.6900
Email: orders@ntis.gov

Cover Photos by Dennis Schroeder: (left to right) NREL 26173, NREL 18302, NREL 19758, NREL 29642, NREL 19795.

NREL prints on paper that contains recycled content.

DISCLAIMER

The Standard Algal Biomass Analytical Methods (Methods) are provided by the National Renewable Energy Laboratory (NREL), which is operated by Alliance for Sustainable Energy, LLC (Alliance) for the U.S. Department of Energy (DOE). These methods were developed and written for commercial research and educational use only.

Access to and use of these Methods shall impose the following obligations on the user. The user is granted the right, without any fee or cost, to use, copy, modify, alter, enhance and distribute these Methods for any purpose whatsoever, except commercial sales, provided that this entire notice appears in all copies of the Methods. The user agrees to credit NREL/Alliance in any publications that result from the use of these Methods. The user also understands that NREL/Alliance is not obligated to provide the user with any support, consulting, training or any training or assistance of any kind with regard to the use of these Methods or to provide the user with any updates, revisions or new versions.

THESE METHODS ARE PROVIDED BY NREL/Alliance "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL NREL/Alliance/DOE BE LIABLE FOR ANY SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES OR ANY DAMAGES WHATSOEVER, INCLUDING BUT NOT LIMITED TO CLAIMS ASSOCIATED WITH THE LOSS OF DATA OR PROFITS, WHICH MAY RESULT FROM AN ACTION IN CONTRACT, NEGLIGENCE OR OTHER TORTIOUS CLAIM THAT ARISES OUT OF OR IN CONNECTION WITH THE ACCESS, USE OR PERFORMANCE OF THESE METHODS.

1. Introduction

- 1.1 This Laboratory Analytical Procedure (LAP) uses a two-step sulfuric acid hydrolysis to hydrolyze the polymeric forms of carbohydrates in algal biomass into monomeric subunits. The monomers are then quantified by either high-performance liquid chromatography (HPLC) or a suitable spectrophotometric method.
- 1.2 An optimized hydrolysis procedure is expected to yield complete hydrolysis of all polymeric structural and storage carbohydrates into monomeric sugars. A range of sequential and optimized inorganic acid hydrolysis conditions with respective hydrolysis liquor collection and analysis should be carried out for algal biomass. Here, only a two-step sulfuric acid hydrolysis and quantification of released monosaccharides are presented, and it is noted that not all carbohydrates may be hydrolyzed in this manner. Future updates to this procedure will involve optimization of the carbohydrate hydrolysis procedures.
- 1.3 Portions of this procedure are substantially similar to ASTM E1758-01 “Standard Method for the Determination of Carbohydrates by HPLC” and the LAP developed for terrestrial feedstocks and reference [1].

2. Scope

- 2.1 This procedure was originally optimized for terrestrial biomass, but has been modified in this LAP to apply to freeze-dried algal biomass.

3. Terminology

- 3.1 *Oven Dry Weight (ODW)*: The weight of the biomass corrected for the percent moisture determined by drying the biomass overnight in a 60°C atmospheric pressure or 40°C vacuum convection oven according to the LAP *Determination of Total Solids and Ash in Algal Biomass* [2].
- 3.2 *Structural Carbohydrates*: The hydrolysable fraction of monomeric sugars after analytical acid hydrolysis with an inorganic acid (such as sulfuric acid), this fraction will also include monosaccharides derived from storage polysaccharides, such as starch. This value should be subtracted from the reporting of total structural carbohydrates.
- 3.3 *Storage Carbohydrates*: The fraction of monomeric glucose measured after enzymatic hydrolysis of algal biomass with an enzyme cocktail of α -amylase and amyloglucosidase. This fraction of storage carbohydrates is referred to as starch, but reflects any portion of α -1,4 linked glucose polymers in algal biomass.
- 3.4 *HPLC*: High-performance liquid chromatography, which separates compounds based on size exclusion and ligand exchange. Details of a comparison between chromatography columns and conditions for HPLC can be found in reference [3].

4. Significance and Use

- 4.1 This procedure was developed for use on terrestrial feedstock and has not been optimized specifically for microalgal biomass. However, the method has been extensively tested and shown to work for algal biomass with no apparent difficulties.
- 4.2 The procedure described here is used to determine the total monomeric carbohydrate content of freeze-dried algal biomass samples after a two-step sulfuric acid hydrolysis.
- 4.3 For the purpose of describing this procedure, we will refer to HPLC quantification, but care should be taken with respect to monosaccharide resolution and chromatography of complex mixtures of algae-specific monosaccharides when choosing an HPLC system. High-performance anion exchange chromatography (HPAEC) may be the method that gives the best separation of major micro-algae-specific monosaccharides [3], but it is not described here.
- 4.4 Alternatively, spectrophotometric quantification of carbohydrates as monosaccharides can be achieved by complexing the free aldehyde group with MBTH (3-methyl-2-benzothiazolinone hydrazone), which has been shown to give an almost equal response of different aldehyde sugars in solution, and matched HPLC quantification of complex mixtures of monosaccharides [4]. This is a novel procedure; a manuscript on the development and validation of this procedure is in preparation and will be included as a reference in future versions of this LAP.
- 4.5 This procedure is used in conjunction with other compositional analysis procedures to determine the summative mass closure for algal biomass.

5. Interferences

- 5.1 Samples with ash content >10% may not be suitable for this procedure as some components of ash may cause side reactions during hydrolysis.
- 5.2 Samples with moisture content >10% may not be suitable for this procedure as the excess moisture will interfere with appropriate acid concentrations.
- 5.3 Samples that are moldy or wet or that have been exposed to an oxygen-rich environment may be compromised, resulting in erroneous carbohydrate values.
- 5.4 This procedure is not suitable for samples containing added acid, base, or catalyst.
- 5.5 Carbohydrates in the hydrolysates should not be measured using a phenol-sulfuric acid quantification procedure. This quantification is known to be susceptible to significant interferences as well as exhibiting substantially different responses for different monosaccharides. Quantification of total carbohydrates by phenol-sulfuric acid failed to match HPLC-quantification of even simple mixtures of sugars.

6. Apparatus

- 6.1 Analytical balance, accurate to 1 mg or 0.1 mg (e.g., Mettler Toledo XP205 DeltaRange)
- 6.2 Vortex mixer
- 6.3 Water bath, set to $30^{\circ}\text{C} \pm 3^{\circ}\text{C}$
- 6.4 Autoclave, suitable for autoclaving liquids, set to $121^{\circ}\text{C} \pm 3^{\circ}\text{C}$
- 6.5 Digital dry block, capable of maintaining $80^{\circ}\text{C} \pm 3^{\circ}\text{C}$ (for spectrophotometric carbohydrate determination) —compatible with glass tubes
- 6.6 Spectrophotometer (for spectrophotometric carbohydrate determination), set to measure at 620 nm
- 6.7 HPLC system equipped with refractive index detector and the following column: Shodex sugar SP0810 (Shodex, #F6378105) with de-ashing cartridge holder (Bio-Rad, #125-0139) and micro-guard de-ashing cartridges (Bio-Rad, #125-0118) or equivalent with corresponding guard column

7. Reagents and Materials Needed

- 7.1 Reagents
 - 7.1.1 Sulfuric acid, 72% w/w (specific gravity 1.6338 at 20°C) (Ricca Chemical Company R8191600-1A)
 - 7.1.2 Calcium carbonate, ACS reagent grade (Sigma Aldrich 239216)
 - 7.1.3 Water, 18.2 MegaOhm ($\text{M}\Omega$) water, 0.2 μm filtered
 - 7.1.4 HPLC method:
 - 7.1.4.1 High purity sugar standards: D(+)glucose (Sigma Aldrich G7528), D(+)xylose (Sigma Aldrich X3877), D(+)galactose (Sigma Aldrich G0750), L(+)arabinose (Sigma Aldrich 10839), D(+)mannose (Sigma Aldrich 63579)
 - 7.1.5 MBTH spectrophotometric method:
 - 7.1.5.1 High purity sugar standards: D(+)glucose (Sigma Aldrich G7528)
 - 7.1.5.2 MBTH (Sigma Aldrich 129739)
 - 7.1.5.3 NaOH, ACS reagent grade (Sigma Aldrich 221465)
 - 7.1.5.4 Dithiothreitol, DTT (Sigma Aldrich G0632)

- 7.1.5.5 Ferric ammonium sulfate dodecahydrate (Fluka 09730)
- 7.1.5.6 Sulfamic acid (Sigma Aldrich 242772)
- 7.1.5.7 Hydrochloric acid (HCl), concentrated (36.5%–38%) (J.T. Baker 9535-03)

7.2 Materials

- 7.2.1 10-mL glass tubes and caps, heavy walled and capable of withstanding high-pressure in an autoclave (VWR, #89079-404). Do not use the caps that are supplied with the vials as they may not allow the samples to vent properly. Rubber stoppers to plug glass tubes (VWR, #15507-144), tear-away aluminum crimp caps (VWR, #16171-851)
- 7.2.2 Adjustable pipet, covering 250 μ L to 10 mL
- 7.2.3 pH paper, range 4–9
- 7.2.4 Cuvettes—glass, quartz, or plastic (spectrophotometric method only)
- 7.2.5 13 x 100-mm glass test tubes or equivalent (spectrophotometric method only)
- 7.2.6 Glass marbles or caps to cover or plug the 13 x 100-mm glass test tubes (7.2.5) (spectrophotometric method only)
- 7.2.7 Disposable syringes (VWR, #BD301073) and compatible 0.2- μ m nylon filters (VWR, #28144-592), if planning to filter as opposed to centrifuging to remove solids
- 7.2.8 1.5-mL capacity HPLC vials (Agilent, #5182-0543) and compatible PTFE/rubber crimp caps (Agilent, #5181-1210) (HPLC method only)
- 7.2.9 Manual or automatic crimpers for caps for glass tubes from 7.2.1 and for HPLC vials from 7.2.8 (if doing the HPLC method)

8. ES&H Considerations and Hazards

- 8.1 Sulfuric acid is corrosive and should be handled with care.
- 8.2 MBTH procedure: hydrochloric acid, sodium hydroxide, and sulfamic acid are moderate health hazards. MBTH and ammonium ferric sulfate dodecahydrate are moderate health hazards. Dithiothreitol is a moderate health hazard with mild flammability and reactivity. Follow all applicable chemical-handling procedures.
- 8.3 Use caution when handling hot pressure tubes after removal from the autoclave. The pressurized tubes are a possible explosion hazard if not cooled properly.

9. Sampling, Test Specimens, and Test Units

- 9.1 The LAP *Determination of Total Solids and Ash in Algal Biomass* [2] should be performed in concurrence with this LAP.
- 9.2 Care must be taken to ensure a representative sample is taken for analysis from an algal biomass prepared and dried by freeze drying, spray drying, etc., ensuring moisture is <10% and is ground/homogenized to a particle size <1mm.

10. Procedure

- 10.1 Prepare the sample for analysis and hydrolyze
 - 10.1.1 Weigh out all samples in concurrence with the LAP *Determination of Total Solids and Ash in Algal Biomass* [2]. The results from this LAP will be used later to correct the biomass weight. The hydrolysis and moisture correction samples must be weighed out on the SAME DAY. If not immediately hydrolyzed, seal samples to prevent any changes in moisture content.
 - 10.1.2 Weigh 25 ± 2.5 mg of freeze-dried algal biomass into a labeled and tared 10-mL glass tube. Record the weight of the sample to the nearest 0.1 mg in a lab notebook.
 - 10.1.3 After all the samples have been weighed out, add 250 μ L of 72% (w/w) sulfuric acid to each tube and vortex to thoroughly mix the acid and biomass. Vortex samples carefully to ensure all solids remain at the bottom of the tube and immersed in the acid.

NOTE: If samples have a tendency to clump, add sulfuric acid to one tube at a time, and vortex the tube immediately after the addition of acid.

NOTE: If the sulfuric acid is tinted yellow or brown, it needs to be replaced.
 - 10.1.4 Place tubes in a $30^{\circ}\text{C} \pm 3^{\circ}\text{C}$ water bath with a water level set to just above the sample/acid level in the tubes. Incubate samples for 1 hour, vortexing each tube vigorously every 5 to 10 minutes. Vortex samples carefully to ensure all solids remain at the bottom of the tube and immersed in acid.
 - 10.1.5 After 1 hour, remove tubes from the water bath and add 7 mL of 18.2 MegaOhm ($M\Omega$) water to each tube. This will bring the sulfuric acid concentration to 4% (w/w).
 - 10.1.6 Place a rubber stopper snugly in the top of the glass tube. Place a tear-away aluminum crimp cap on top of the rubber stopper. Crimp the caps tightly and vortex to mix the contents thoroughly. Tubes may also be inverted to mix, ensure that the solids remain in the liquid, and do not end up stuck to the rubber stopper.

- 10.1.7 Place the tubes in a rack suitable for autoclaving. Be sure to place the rack in a tray for secondary containment. Secondary containment should contain any tubes and their contents if they were to break during the autoclave cycle.
- 10.1.8 Autoclave the sealed samples for one hour at 121°C using the liquids (slow) setting.
- 10.1.9 After completion of the autoclave cycle, open the autoclave door, but allow tubes to cool for approximately 15 minutes before removing them from the autoclave.
- 10.1.10 Once the samples are removed from the autoclave, allow the samples to cool for 30 minutes to 1 hour or until tubes reach room temperature.

10.2 Prepare samples for HPLC analysis

- 10.2.1 Hydrolyzed samples for HPLC analysis should be neutralized the same day or the solids separated from the hydrolysate and the hydrolysate (still acidic) may be stored in a refrigerator (4°C) for up to 24 hours before being neutralized. Watch for precipitate formation after the sample has been refrigerated. The sample must be vortexed thoroughly, and then re-filtered or centrifuged to remove the precipitate before analysis by HPLC.

NOTE: If longer storage is necessary, store the hydrolysate (separated from the solids) in a freezer (-20°C) for up to 7 days. DO NOT neutralize before freezing.

- 10.2.2 Once the samples are cooled to room temperature, vortex each sample and take an aliquot for neutralization. Aliquots are typically about 1–3 mL; do not attempt to neutralize large volumes of hydrolysate at once. The sample does not need to be filtered or centrifuged to remove the solids before aliquoting.
- 10.2.3 Place the aliquot in a vessel suitable for neutralizing the sample. The vessel must be large enough to contain the sample when it off-gases and bubbles during the neutralization process. Examples might be 25–50-mL Erlenmeyer flasks or 50-mL plastic centrifuge tubes (preferable to remove solids and precipitate from hydrolysate if centrifuging).
- 10.2.4 Neutralize the aliquot to a pH between 6 and 8 using calcium carbonate. It is very important to neutralize the sample correctly as low and high pHs can ruin an HPLC column. The amount of calcium carbonate required to neutralize the sample will depend on several factors, including, but not limited to, aliquot amount and sample buffering capacity (e.g., ash content). **Add calcium carbonate in small amounts** and vortex between additions. Allow the sample sufficient time to finish bubbling before checking pH and/or adding any additional calcium carbonate. Check the pH periodically during the addition of the calcium carbonate. DO NOT add large amounts of

calcium carbonate at one time, this often leads the sample to froth and turn into a lump of precipitate. A precipitate will form as the neutralization process progresses; this is normal.

10.2.5 The neutralized sample should be filtered through a 0.2- μm nylon filter (use the disposable syringes and compatible filters) OR centrifuged to remove ALL solids and precipitate. For the HPLC method, place the filtered or centrifuged hydrolysate in a clean, labeled HPLC vial (1.5-mL vials) and cap tightly with a crimp cap (caps are considered tight when they cannot be turned by hand).

10.3 Prepare samples for spectrophotometric analysis

10.3.1 Hydrolyzed samples for spectrophotometric analysis should be prepared by separating the solids from the hydrolysate. Acidic solids-free hydrolysate may be stored in a refrigerator (4°C) for up to 24 hours before being analyzed.

NOTE: If longer storage is necessary, store the hydrolysate (separated from the solids) in a freezer (-20°C) for up to 7 days. DO NOT neutralize before freezing.

10.3.2 The sample should be filtered through a 0.2- μm nylon filter (use the disposable syringes and compatible filters) OR centrifuged to remove ALL solids. For the spectrophotometric method, place the filtered or centrifuged hydrolysate in a clean, labeled vial and cap.

10.4 HPLC analysis for monomeric sugars

10.4.1 Prepare a series of calibration standards containing the compounds to be quantified. Refer to Tables 1 and 2 for suggested concentration ranges and sugars. Use at least a four-point calibration, for example, 0.05, 0.25, 1.25, and 4 mg/mL with a 2.5 mg/mL calibration verification standard (CVS) according to the dilution table below:

Table 1. Suggested Dilutions for Calibration Curve Sample Preparations

Concentration of stock (mg/mL)	Dilution (stock + water)	Final concentration (mg/mL)
10	400 + 600	4
10	125 + 875	1.25
1	250 + 750	0.25
1	50 + 950	0.05
10	250 + 750	~2.5 (CVS)

Table 2. Suggested Standard Concentrations for the HPLC. Choice of monosaccharides is based on SP0810 column chromatography and is not an accurate reflection of algal monosaccharides. (For a full list of algae-specific monosaccharides and suggested calibration mixtures, see reference [3]).

Component	Suggested concentration range (mg/mL)
D(+)-glucose	0.05–4
D(+)-xylose	0.05–4
D(+)-galactose	0.05–4
L(+)-arabinose	0.05–4
D(+)-mannose	0.05–4
CVS	~2.5

10.4.2 Prepare an independent CVS for each calibration set. The CVS should come from a different sugar stock solution than the standards. The concentration of the CVS should fall somewhere in the middle of the chosen calibration range. The CVS should be run immediately after the calibration, after every 10 samples, and at the end of the HPLC run. The CVS is used to verify the quality and stability of the calibration curve throughout the run.

10.4.3 Analyze the calibration standards, CVS, and samples by HPLC using a Shodex sugar SP0810 column equipped with an appropriate guard column (keep the guard column out of the column compartment; it should be kept at room temperature) and a refractive index detector. Be sure to flush both the column and guard column if they are new. Flush them separately into a waste container. DO NOT flush them into the detector or flush the guard column into the sugar separation column. Running at the method temperature and maximum flow rate for approximately 1 hour should be sufficient:

- Sample volume: 50 μ L
- Mobile phase: 18.2 MegaOhm (M Ω) water, 0.2 μ m filtered
- Flow rate: 0.6 mL/min
- Column temperature: 85°C
- Detector temperature: 55°C
- Run time: 50 minutes

NOTE: Run time may be extended if there are late eluting compounds and/or the addition of post-run time may be necessary to elute sugar breakdown products from the column.

10.4.4 Quantify the monomeric sugars by Chemstation (Agilent) or equivalent LC software. Correlation for each calibrated compound should be 0.9999 or better.

10.5 Spectrophotometric analysis for monomeric sugars using MBTH

10.5.1 Preparation of reagents/solutions—use 18.2 MegaOhm (M Ω) water for all preparations requiring water:

- 0.5M NaOH solution—store at room temperature in a plastic (HDPE) container and remake once a precipitate forms (looks like silvery flakes) — typically about 2–3 weeks
- 3 mg/mL MBTH—store in the dark at 4°C for up to one week
- 1mg/mL DTT—store in the dark at 4°C for up to one week
- 0.25M HCl—prepare in a 200-mL volumetric flask and store at room temperature in the dark for up to one month:
 - Add a small amount of 18.2 M Ω water to the volumetric flask
 - Pipet 4.16 mL of concentrated HCl (~36.5%-38%) into the flask
 - Bring to volume with 18.2 M Ω water
 - Cap and mix thoroughly
- Ferric solution: 0.5% ferric ammonium sulfate dodecahydrate, 0.5% sulfamic acid (w/v), and 0.25M HCl. Dissolve 200 mg of ferric ammonium sulfate dodecahydrate and 200 mg of sulfamic acid into 40 mL of 0.25M HCl. This solution may be stored at room temperature in the dark for up to one month.

10.5.2 Preparation of working solution, made up fresh right before adding to samples:

- MBTH working solution: 1:1 (v/v) of 3 mg/mL MBTH and 1 mg/mL DTT

10.5.3 Preparation of standards:

- Prepare glucose standard and CVS stock solutions (these must be made up separately):
 - Weigh approximately 25 mg of glucose into a clean, tared, class A, 100-mL volumetric flask
 - Record the weight to the nearest 0.1 mg

- Bring to volume with 18.2 MΩ water
- Cap and mix thoroughly
- The standard and CVS working solution concentrations should be approximately 0.25 mg/mL glucose. Due to balance limits, start with a higher stock glucose concentration and dilute to 0.25 mg/mL. ALWAYS prepare the standards and CVS from separate stocks so that you can check the accuracy of your calibration with the CVS.
- Glucose stock solutions should be stored in a freezer (-20°C) for up to 2 months or 5 days in the refrigerator (4°C).
- Create a six-point calibration with the 0.25 mg/mL glucose stock solution using the dilutions outlined in Table 3. Prepare dilutions directly into the glass reaction vials (13 x 100 mm or equivalent).

Table 3. Suggested Standard Concentrations for the Spectrophotometric Method

Concentration of stock/working (mg/mL)	Dilution stock + water (μL)	Final concentration (mg/mL)
0.25	0 + 500	0
0.25	20 + 480	0.010
0.25	30 + 470	0.015
0.25	50 + 450	0.025
0.25	75 + 425	0.0375
0.25	100 + 400	0.050

- Prepare the CVS at the center of the calibration range: 75 μL CVS glucose stock solution (0.25 mg/mL) and 425 μL 18.2 MegaOhm (MΩ) water. Prepare the CVS dilution directly into the glass reaction vials.

10.5.4 Reaction

- Preheat a digital dry block that fits 13 x 100-mm glass test tubes (or equivalent reaction vials) to 80°C ± 3°C.

Prepare sample dilutions directly in the glass reaction vials: start with a 1:10 (50 μL hydrolysate) or greater dilution of the filtered/centrifuged acidic hydrolysate. Dilute the sample to a total volume of 0.5 mL. If the sample is known to be high in carbohydrates, start with a 1:20 dilution. The goal is to have the sample absorbance fall in the middle of the calibration curve.

NOTE: DO NOT use neutralized hydrolysate as the soluble salts leftover from the neutralization will affect the reaction.

NOTE: DO NOT use less than a 1:10 dilution on acidic hydrolysates as the excess acid can have an effect on the reaction.

- ALWAYS run a set of standards and a CVS with each batch of samples.
- Add 500 μL 0.5M NaOH to each glass tube.
- Add 500 μL MBTH working solution (section 10.5.2) to each glass tube, vortex carefully to mix, and cover the tubes with a glass marble or cap. Immediately place the tubes with solutions in the preheated dry block at 80°C for 15 ± 1 min. Do not exceed the incubation time. Work in smaller sample sets (<20). Do not allow the samples and standards to sit around after addition of the NaOH or MBTH working solution—they must be placed on the hot block within a few minutes of adding these solutions.
- After the 15-minute incubation time, turn off the block and immediately add 1 mL of the ferric solution while the glass tubes are still on the block (section 10.5.1).
- Once the ferric solution has been added, remove the glass tubes from the hot block and carefully vortex to mix.
- Allow the samples and standards to react with the ferric solution for 10–15 minutes while the samples cool to room temperature. Tubes may remain open during this time. Allow them to cool in a hood to control fumes.
- Once the samples are at room temperature, add 2.5 mL 18.2 M Ω water, mix by pipetting or vortexing.
- Once mixed, place an aliquot of the samples and standards into the appropriate cuvettes and obtain an absorbance on the spectrophotometer at 620 nm. Take an absorbance within an hour after the final dilution step. Zero the spectrophotometer on 18.2 M Ω water. DO NOT zero the spectrophotometer on the 0 standard from the calibration standard set.
- Use the glucose calibration curve and linear regression to quantify the total carbohydrate concentration in mg/mL. Be sure to include the 0 standard as a point in the calibration. Remember to correct the sample carbohydrate concentration for the amount of hydrolysate that was used in the dilution (for a 1:10, use 50 μL) and then get a total carbohydrate content by multiplying that concentration by the total volume from the hydrolysis—7.25 mL. When building a linear calibration curve, DO NOT force the curve through the origin.

11. Calculations

- 11.1 Calculate the ODW of the sample, using the average total solids content as determined by the LAP *Determination of Total Solids and Ash in Algal Biomass* [2]:

$$ODW = \frac{Weight_{sample} \times \% Total Solids}{100}$$

- 11.2 Calculate monomeric sugar concentrations (mg/mL) for each sample using the linear regression coefficients for the spectrophotometric glucose determination or the HPLC calibration.
- 11.3 Calculate monomeric sugar content (mg) for each sample using the following equation:

$$Total\ Monosaccharides_{(mg)} = Monosaccharide_{(mg/mL)} \times 7.25_{(mL)}$$

- 11.4 Sum the monomeric sugars to get total sugars in mg (for the HPLC method).
- 11.5 Calculate the amount of monomeric (structural) sugar in the sample on a percent ODW basis:

$$\% Monomeric\ Sugar = \left(\frac{Total\ Monosaccharides_{(mg)}}{ODW_{Sample}} \right) - \% Starch \times 100$$

- 11.6 To report or calculate the relative percent difference (RPD) between two samples, use the following calculation:

$$RPD = \left(\frac{(X_1 - X_2)}{X_{mean}} \right) \times 100$$

where:

X_1 and X_2 = measured values

X_{mean} = the mean of X_1 and X_2

- 11.7 To report or calculate the root mean square deviation (RMS) or the standard deviation (STDEV) of the samples, use the following calculation:

$$RMS = x_m = mean = \sqrt{\left(\frac{\sum x}{n} \right)^2}$$

$$RMS\ deviation = \sigma = stdev = \sqrt{\frac{\sum_{i=1}^n (x_i - x_m)^2}{n}}$$

where:

x_m = the root mean square of all x values in the set

n = number of samples in set

x_i = measured value from the set

12. Report Format

- 12.1 Report the total carbohydrate results as a percent of the ODW. Individual sugars may also be reported as a percent of the ODW.
- 12.2 For replicate analyses of the same sample, report the average, standard deviation, and %RPD.
- 12.3 Report the CVS recovery.

13. Precision and Bias

- 13.1 Precision and bias need to be determined by a round robin experiment using standard method verification biomass. The target RPD should be less than 10%.

14. Quality Control

- 14.1 *Reported Results*: Report results with two decimal places. Report the average, standard deviation, and %RPD.
- 14.2 *Replicates*: Run all samples in triplicate, unless prohibited by the amount of sample available (hydrolysis only).
- 14.3 *RPD Criterion*: Determined by data quality objectives and laboratory-specific *Quality Assurance Plan*.
- 14.4 *CVS*: CVSs should be independently prepared and analyzed as per the procedure. Required agreement for calibration verification standard quantification relative to the theoretical concentration should be within 5% RPD.
- 14.5 *Sample Size*: 25 ± 2.5 mg for each replicate.
- 14.6 *Sample Storage*: Hydrolysis liquors may be separated from the acid insoluble solids and stored (still acidic) for up to two weeks in a refrigerator (4°C). Neutralized liquors may be stored in a refrigerator (4°C) for up to 4 days.
- 14.7 *Standard Storage*: HPLC standards should be stored in a freezer (-20°C) and removed when needed. Thaw and vortex standards prior to use. Filter if necessary.

14.8 *Standard Preparation*: Standards should be prepared as described in the procedure.

14.9 *Definition of a Batch*: Any number of samples that are analyzed and recorded together.

14.10 *Control Charts*: Quality assurance/quality control material should be control charted to verify reproducibility.

15. Appendices

15.1 List of revisions/updates:

- Distribution of May 16, 2013 DRAFT version
- Revision July 8, 2013, addition of additional experimental details and dilution descriptions
- Revision July 17, 2013, correction to dilutions in Table 2
- Revision July 26, 2013, addition of CVS recovery target precision
- Revision December 29, 2015, update of spectrophotometric and HPLC methods and sample storage

16. References

[1] ASTM E1758-01 “Standard Method for the Determination of Carbohydrates by HPLC.” *2003 Annual Book of ASTM Standards*, Volume 11.05. Philadelphia, PA: American Society for Testing and Materials, International.

[2] Van Wychen, S.; Laurens, L.M.L. *Determination of Total Solids and Ash in Algal Biomass*. NREL/TP-5100-60956. Golden, CO: National Renewable Energy Laboratory. 2015. <http://www.nrel.gov/docs/fy14osti/60956.pdf>.

[3] Templeton, D.W.; Quinn, M.; Van Wychen, S.; Hyman, D.; Laurens, L.M.L. “Separation and quantification of microalgal carbohydrates.” *J. Chrom. A* (1270), 2012; pp. 225–234.

[4] Anthon, G.E.; Barrett, D.M. “Separation and quantification of microalgal carbohydrates.” *Anal. Biochem.* (305), 2002; pp. 287–289.