



# Determination of Carbonyls in Pyrolysis Bio-oils by Potentiometric Titration: Faix Method

## Laboratory Analytical Procedure (LAP)

**Issue Date: March 2, 2016**

Stuart Black and Jack Ferrell  
*National Renewable Energy Laboratory*

Mariefel V. Olarte and Asanga B. Padmaperuma  
*Pacific Northwest National Laboratory*

**NREL is a national laboratory of the U.S. Department of Energy  
Office of Energy Efficiency & Renewable Energy  
Operated by the Alliance for Sustainable Energy, LLC**

This report is available at no cost from the National Renewable Energy  
Laboratory (NREL) at [www.nrel.gov/publications](http://www.nrel.gov/publications).

**Technical Report**  
NREL/TP-5100-65888  
March 2016

Contract No. DE-AC36-08GO28308

# Determination of Carbonyls in Pyrolysis Bio-oils by Potentiometric Titration: Faix Method

## Laboratory Analytical Procedure (LAP)

**Issue Date: March 2, 2016**

Stuart Black and Jack Ferrell  
*National Renewable Energy Laboratory*

Mariefel V. Olarte and Asanga B.  
Padmaperuma  
*Pacific Northwest National Laboratory*

Prepared under Task Number BB143810

**NREL is a national laboratory of the U.S. Department of Energy  
Office of Energy Efficiency & Renewable Energy  
Operated by the Alliance for Sustainable Energy, LLC**

This report is available at no cost from the National Renewable Energy  
Laboratory (NREL) at [www.nrel.gov/publications](http://www.nrel.gov/publications).

## NOTICE

This report was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or any agency thereof.

This report is available at no cost from the National Renewable Energy Laboratory (NREL) at [www.nrel.gov/publications](http://www.nrel.gov/publications).

Available electronically at SciTech Connect <http://www.osti.gov/scitech>

Available for a processing fee to U.S. Department of Energy and its contractors, in paper, from:

U.S. Department of Energy  
Office of Scientific and Technical Information  
P.O. Box 62  
Oak Ridge, TN 37831-0062  
OSTI <http://www.osti.gov>  
Phone: 865.576.8401  
Fax: 865.576.5728  
Email: [reports@osti.gov](mailto:reports@osti.gov)

Available for sale to the public, in paper, from:

U.S. Department of Commerce  
National Technical Information Service  
5301 Shawnee Road  
Alexandria, VA 22312  
NTIS <http://www.ntis.gov>  
Phone: 800.553.6847 or 703.605.6000  
Fax: 703.605.6900  
Email: [orders@ntis.gov](mailto:orders@ntis.gov)

*Cover Photos by Dennis Schroeder: (left to right) NREL 26173, NREL 18302, NREL 19758, NREL 29642, NREL 19795.*

NREL prints on paper that contains recycled content.

## **DISCLAIMER**

The Standard Bio-oil Analytical Methods (Methods) are provided by the National Renewable Energy Laboratory (NREL), which is operated by Alliance for Sustainable Energy, LLC (Alliance) for the U.S. Department of Energy (DOE). These methods were developed and written for commercial research and educational use only.


Access to and use of these Methods shall impose the following obligations on the user. The user is granted the right, without any fee or cost, to use, copy, modify, alter, enhance and distribute these Methods for any purpose whatsoever, except commercial sales, provided that this entire notice appears in all copies of the Methods. The user agrees to credit NREL/Alliance in any publications that result from the use of these Methods. The user also understands that NREL/Alliance is not obligated to provide the user with any support, consulting, training or any training or assistance of any kind with regard to the use of these Methods or to provide the user with any updates, revisions or new versions.

THESE METHODS ARE PROVIDED BY NREL/Alliance "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL NREL/Alliance/DOE BE LIABLE FOR ANY SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES OR ANY DAMAGES WHATSOEVER, INCLUDING BUT NOT LIMITED TO CLAIMS ASSOCIATED WITH THE LOSS OF DATA OR PROFITS, WHICH MAY RESULT FROM AN ACTION IN CONTRACT, NEGLIGENCE OR OTHER TORTIOUS CLAIM THAT ARISES OUT OF OR IN CONNECTION WITH THE ACCESS, USE OR PERFORMANCE OF THESE METHODS.

## 1. Introduction

- 1.1 While pyrolysis bio-oils are comprised of a large variety of compounds and chemical functional groups, quantification of carbonyl groups is especially important. Carbonyls are known to be responsible for the instability of bio-oil during both storage and processing. The titration method presented here is a simple technique which can reliably quantify the total carbonyl content of bio-oils.
- 1.2 For analysis of bio-oils, quantification of carbonyl groups by titration has traditionally been accomplished using the method of Nicolaides [1]. This is a simple procedure where carbonyls are converted to the corresponding oxime (see Scheme 1). The liberated HCl reacts with pyridine to force the equilibrium to completion. The conjugate acid of pyridine is titrated with a known amount of NaOH (base titrant). The number of equivalents of NaOH used is stoichiometrically equivalent to the moles of carbonyl present in the bio-oil.

Scheme 1:



- 1.3 The Nicolaides method, however, has several limitations. It can require reaction times in excess of 48 hours to reach completion. This severely limits sample throughput. It utilizes pyridine, which is toxic. Sample weights of 1 to 2 g are required. Sample weight used is dependent on the amount of hydroxylamine HCl present and the carbonyl content of the sample. If initial estimates of the sample weight used are incorrect, the titration has to be repeated. Furthermore, it has recently been shown that the Nicolaides method significantly underestimates carbonyl content of bio-oils [3].
- 1.4 Faix, et al. [2] developed a method that has been modified here to address the issues of the Nicolaides method. The reaction is carried out at 80 °C for 2 hours, thereby increasing sample throughput. Pyridine has been replaced with triethanolamine which is a less toxic chemical. The sample size can be reduced to 100 to 150 mg. The triethanolamine consumes the liberated HCl, driving the reaction to completion and the unconsumed triethanolamine is titrated directly. A secondary titration of the hydroxylamine is unnecessary. An in-depth comparison of both of these carbonyl titration methods has recently been published [3].

## 2. Scope

- 2.1 The method described here has been modified from the original method [2] to be more applicable to the analysis of pyrolysis bio-oils. This method was developed for the analysis of raw pyrolysis bio-oils, but it has been successfully applied to other types of

biomass-derived oils, including hydrotreated bio-oils. Additionally, this method has been used to monitor changes in carbonyl content during both aging and upgrading.

### 3. Terminology

- 3.1 *Bio-oil* – The crude liquid product of converting solid biomass into a liquid via fast pyrolysis or other thermochemical conversion process.
- 3.2 *Carbonyl* – the chemical functional group consisting of a carbon-oxygen double bond, C=O. For this method, this includes all aldehydes and ketones, but does not include carboxylic or lactone groups.

### 4. Interferences

- 4.1 The selectivity of the method was tested by using 1-butanol, 1-pentanol, tertiary butanol, 2-propanol, ethyl acetate, acetic acid, xylose and glucose as model compounds, representing alcohol, ester, carboxylic acid and carbohydrates in the bio-oil.
- 4.2 No interferences were seen for ethyl acetate or acetic acid. Monosaccharides are measured using this method. Addition of alcohols cause interferences but it is dependent on chain length. The reason is as yet undetermined but may be related to solvent properties of the alcohol.

### 5. Apparatus

- 5.1 Analytical balance accurate to 0.1 mg
- 5.2 Dry block heater with magnetic stirrer or hot water bath with magnetic stirrer
- 5.3 Automatic titration instrument (Metrohm Titrando 809 or equivalent)

### 6. Reagents and Materials Needed

- 6.1 Reagents
  - 6.1.1 Deionized water
  - 6.1.2 Ethanol (reagent grade) (CAS # 64-17-5)
  - 6.1.3 Hydroxylamine hydrochloride (CAS # 5470-11-1)
  - 6.1.4 Triethanolamine (CAS #102-71-6)
  - 6.1.5 Hydrochloric acid (37%) (CAS # 7647-01-0)
  - 6.1.6 Sodium Carbonate (primary standard) (CAS # 497-19-8)
  - 6.1.7 4-(benzyloxy)benzaldehyde (CAS # 4397-53-9)

- 6.1.8 Dimethyl sulfoxide (CAS # 67-68-5)
- 6.2 Materials
  - 6.2.1 5 mL glass Reacti-vials (Thermoscientific TS-13223) with solid lid
  - 6.2.2 Teflon spinvane for 5 ml Reacti-vial
  - 6.2.3 200 mL volumetric flask
  - 6.2.4 Volumetric or mechanical pipettes
- 6.3 Reagent solutions
  - 6.3.1 *Hydroxylamine hydrochloride solution (Solution A)*: Add 7.7 g of hydroxylamine hydrochloride and 50 mL of deionized water to a 250 mL volumetric flask. Swirl till all solids are dissolved, then dilute up to the mark with ethanol.
  - 6.3.2 *Triethanolamine solution (Solution B)*: Add 17.4 mL of triethanolamine to a 250 mL volumetric flask, then dilute up to the mark with ethanol.
  - 6.3.3 *Hydrochloric acid solution*: Either purchase 0.1N solution or prepare using 10 ml concentrated HCl and 1 L water.

## 7. ES&H Considerations and Hazards

- 7.1 Ethyl alcohol is flammable.
- 7.2 Follow all applicable chemical handling procedures.
- 7.3 Follow all applicable waste disposable and handling procedures.

## 8. Sampling, Test Specimens and Test Units

- 8.1 Make sure the oil sample is at room temperature prior to withdrawing a sample. Bio-oil should be thoroughly homogenized to obtain a representative sample.
- 8.2 The test specimen should consist of 100 to 150 milligrams of bio-oil.
- 8.3 Exposure to oxygen and heat should be minimized to prevent sample degradation prior to analysis.

## 9. Analytical Procedure

- 9.1 Standardization of the base solution.
  - 9.1.1 Prepare primary standard sodium carbonate using standard methods.

- 9.1.2 Weigh 100 to 150mg of sodium carbonate to a titration vessel, record the actual weight, add a stir bar and add appropriate amount of D.I. water.
  - 9.1.3 Quantitatively transfer the sample to a titration vessel and titrate with the acid solution using an automatic titrator, and record the endpoint.
  - 9.1.4 Repeat the process twice, to obtain three points.
  - 9.1.5 Use the average value as the normality of the acid solution.
- 9.2 Preparation of titration blanks
- 9.2.1 Blank A: add 1 mL dimethyl sulfoxide (DMSO) to 5 ml Reacti-vial with spinvane.
  - 9.2.2 Add 2 mL hydroxylamine hydrochloride solution (solution A).
  - 9.2.3 Add 2 ml triethanolamine solution (solution B).
  - 9.2.4 Cap tightly, place in preheated (80 °C) heater block or water bath and stir for 2 hours.
  - 9.2.5 Titrate with acid solution using an automatic titrator and record endpoint.
  - 9.2.6 Blank B: If mineral acid is suspected to be present in sample, add 1 mL DMSO to 5 mL Reacti-vial with spin vane.
  - 9.2.7 Add 2 mL triethanolamine solution (solution B).
  - 9.2.8 Cap tightly and stir at 80 °C for 2 hours.
  - 9.2.9 Quantitatively transfer the sample to a titration vessel and titrate with acid solution and record endpoint.
  - 9.2.10 Repeat process three times to obtain three points.
- 9.3 Validation of the method using a known carbonyl.
- 9.3.1 Weigh ~ 100mg of 4-(benzyloxy)benzaldehyde (4-BBA) to a Reacti-vial, record the actual weight, add a spinvane.
  - 9.3.2 Add 1 mL DMSO.
  - 9.3.3 Dissolve the sample in 2 mL hydroxylamine hydrochloride solution (solution A).
  - 9.3.4 Add 2 ml triethanolamine solution (solution B).



- 9.3.5 Close lid tightly and stir at 80 °C for 2 hours.
  - 9.3.6 Quantitatively transfer the sample to a titration vessel with ethanol and water to make a solution of 80% ethanol and titrate with the acid solution using an automatic titrator, and record the endpoint.
  - 9.3.7 Repeat the process three times, to obtain three points.
- 9.4 Analysis of bio-oil using the method
- 9.4.1 Weigh close to 100 mg of the bio-oil to Reacti-vial, record the actual weight and add a spinvane.
  - 9.4.2 Add 1 mL DMSO.
  - 9.4.3 Dissolve the sample in 2 mL hydroxylamine hydrochloride solution (solution A).
  - 9.4.4 Add 2 ml triethanolamine solution (solution B).
  - 9.4.5 Close lid tightly and stir at 80 °C for 2 hours.
  - 9.4.6 Quantitatively transfer the sample to a titration vessel with ethanol and water to make a solution of 80% ethanol and titrate with the acid solution using an automatic titrator, and record the endpoint.
  - 9.4.7 Repeat the process three times, to obtain three points.
  - 9.4.8 Blank C: If mineral acid is suspected to be present in sample, Weigh close to 100 mg of the bio-oil into a Reacti-vial, record weight and add a spinvane.
  - 9.4.9 Add 1 mL DMSO.
  - 9.4.10 Add 2 mL triethanolamine solution (solution B).
  - 9.4.11 Cap tightly and stir at 80 °C for 2 hours.
  - 9.4.12 Quantitatively transfer the sample to a titration vessel with ethanol and water to make a solution of 80% ethanol and titrate with acid solution and record endpoint.
  - 9.4.13 Repeat process three times to obtain three points.
- 9.5 Other considerations
- 9.5.1 Blanks B and C are unnecessary if only organic acids are present in the bio-oil sample.

- 9.5.2 4-(benzyloxy)benzaldehyde (4-BBA) has a sufficiently high molecular weight to use approximately the same amount as a bio-oil sample. Other models may be used.
- 9.5.3 An 80% solution of ethanol is not necessary for the transfer of the sample to the titration vessel. Only a final volume of 80% ethanol is needed.
- 9.5.4 Titration of the sample following oximation should be done within 2 to 8 hours. Triethanolamine can form triethanolamine•HCl which will result in an error in the measurement.
- 9.5.5 Faix, et al. [2] mix the triethanolamine and hydroxylamine hydrochloride solutions to make a stock solution. This leads to the formation of triethanolamine•HCl which will result in a measurement error.

## 10. Results

10.1 The following Tables can be used as a guide to record the data.

10.2 Standardize the base solution:

Purity of Na <sub>2</sub> CO <sub>3</sub>		%
	Na <sub>2</sub> CO <sub>3</sub> weight (g)	Endpoint (mL)
Run #1		
Run #2		
Run #3		

10.3 Blank endpoint volume. Repeat as necessary for all blanks:

	Endpoint (mL)
Run #1	
Run #2	
Run #3	

10.4 Validation of the method using 4-(benzyloxy)benzaldehyde (4-BBA):

Purity of 4-BBA		%
	4-BBA (g)	EP(mL)
Run #1		
Run #2		
Run #3		

## 10.5 Analysis of bio-oils:

	Bio-oil weight (g)	Endpoint (mL)
Run #1		
Run #2		
Run #3		

10.5.1 A typical titration curve consists of single endpoint.

10.5.2 If the software is incapable of detecting the endpoint, the 1st derivative of the titration curve can be constructed using a standard graphing program.

## 11. Calculations

11.1 Standardization of the base solution.

11.1.1 If the weight of the dry sodium carbonate in grams is  $w_1$ , the purity (written as a fraction, i.e., 99% is 0.99), and the endpoint in mL, the concentration of the acid solution (mol/L) is given by:

$$[acid] = (1000 * w_1 * \text{purity of Na}_2\text{CO}_3) / (105.9885 * \text{endpoint})$$

11.2 Validation of the method using a known carbonyl.

11.2.1 If the weight of the 4BBA in grams is  $w_2$ , the purity (written as a fraction, i.e., 99% is 0.99) is known, the concentration of the acid solution is [acid], triethanolamine/hydroxylamine•HCl blank is  $EP_{BA}$  (the average value of three blanks, in mL) and the endpoint is EP (in mL); the concentration of 4BBA (mol/L) in the sample is given by:

$$[4BBA] = \left( \frac{(EP_{BA} - EP)}{w_2} \right) * [acid]$$

11.3 Analysis of bio-oil

11.3.1 If the weight of the bio-oil in grams is  $w_3$ , the concentration of the acid solution is [acid], triethanolamine/hydroxylamine•HCl blank is  $EP_{BA}$  (the average value of three blanks in mL) and the endpoint is EP (in mL); the concentration of carbonyls in bio-oils [CO] (mmol/g-bio oil) is given by:

$$[CO] = \left( \frac{(EP_{BA} - EP)}{w_3} \right) * [acid]$$

11.4 Acid correction

- 11.4.1 The presence of mineral acids or organic acids with  $pK_a < 2$  can cause artificially low carbonyl values due to the reaction of the acid with triethanolamine. If this is suspected, blanks detailed in sections 9.2.6 and 9.4.8 should be performed. The weight of bio-oil in the sample in grams is  $w_3$ , EP is the endpoint of the sample and  $EP_{BA}$  is the triethanolamine/hydroxylamine blank.  $EP_{BB}$  is the endpoint of Blank B,  $EP_{BC}$  is endpoint of Blank C and the weight (g) of oil used in Blank C is  $w_{BC}$ :

$$[CO_{corrected}] = \left( \frac{(EP_{BA} - EP)}{w_3} - \frac{EP_{BB} - EP_{BC}}{w_{BC}} \right) * [acid]$$

For bio-oil samples that contain mostly acetic acid, this is an unnecessary step.

## 12. Report Format

- 12.1 Report the concentration of carbonyls in bio-oils as mmol of CO per gram of bio-oil. The standard deviation may also be reported.

## 13. Precision and Bias

- 13.1 In 2015, an inter-laboratory study was performed on a raw pyrolysis bio-oil using the method as described here [4]. It was found that the Faix carbonyl method presented here was especially reliable, with < 5% RSD variability between all participating laboratories.

## 14. Quality Control

- 14.1 To ensure proper operation, validation with a known carbonyl (9.3) should be performed periodically.
- 14.2 Replicates: Run all samples in triplicate.

## 15. References

[1] Nicolaides, G. M. (1984). *The chemical characterization of pyrolytic oils*. MASc Thesis, Dept. of Chemical Engineering, University of Waterloo, Waterloo, Ontario, Canada, (1984)

[2] Faix, O., B. Andersons and G. Zakis, *Determination of Carbonyl Groups of Six Round Robin Lignins*. *Holzforschung*, 1998. 52: p. 268-272

[3] Black, S. and Ferrell, J. *Determination of Carbonyl Groups in Pyrolysis Bio-oils Using Potentiometric Titration: Review and Comparison of Methods*. *Energy & Fuels*, 2016. DOI: 10.1021/acs.energyfuels.5b02511

[4] Ferrell, J.R. III; Olarte, M.V; Christensen, E.D.; Padmaperuma, A.B.; Connatser, R.M.; Stankovikj, F.; Meier, D.; Paasikallio, V. *Standardization of Chemical Analytical Techniques for Pyrolysis Bio-oil: History, Challenges, and Current Status of Methods*. Biofuels, Bioproducts & Biorefining 2015, submitted.